

Tanita Body Composition Analyser

Model : TIMC980

Body composition analyser featuring state of the art technology including multi frequency analysis and colour touch screen. Ideal for all medical practitioners, hospitals and health clubs.

FEATURES

TIMC980

- state of the art technology including multi frequency analysis, colour touch screen technology and Microsoft® Windows® real time operating system
- segmental analysis for each leg, arm and trunk including body balance, fat distribution and reactance/resistance assessment allowing in-depth trend consultation
- 2 user modes: standard (5–99 years of age) athlete (18–99 years of age), both male and female with height range of 90–250 cm
- user friendly and designed to be used as a stand alone unit allowing clients to take a measurement without assistance thanks to the extra large interactive touch screen display
- six frequencies allow heightened accuracy using the latest Tanita equations
- increased maximum weighing capacity of 300 kg
- full segmental body composition analysis taken in under 30 seconds
- a modular system for convenient transportation. The interlocking system can be set up in under 5 minutes
- Australian Register of Therapeutic Goods (ARTG) #142294

SPECIFICATIONS

Capacity:	300 kg x 100 g
Platform size (mm):	450 (W) x 490 (D) x 65 (H)
Net Weight:	33 kg

OPTIONS

- GMON health management software, charts patient's progress



Specifications may change without notice