

# Portable Height Measuring Rod / Stadiometer

Model : WMHM200P

The WMHM200P mechanical stadiometer disassembles into small pieces for easy transportation using the included carry case. Cleverly designed to enable height measurement of people standing up or lying down which makes it suitable for measuring baby height as well as adults.

It can be set up nearly anywhere for measurement at home, in offices, or in clinics.


## FEATURES

WMHM200P

- stabilisers to help stabilise the stadiometer when measuring against a wall
- portable for professionals on the move
- can be adjusted to measure height of patients in a lying down position
- disassembles to fit in the supplied carry case
- Australian Register of Therapeutic Goods (ARTG) # 329078

## SPECIFICATIONS

Measuring range:	14 cm - 205 cm
Graduation:	1 mm
Overall dimensions (mm):	350 (W) x 400 (D) x 2430 (H)
Net weight:	1.9 kg



Horizontal measurement mode  
(patient lying down)



Carry case



Vertical measurement mode  
(patient standing)

Specifications may change without notice