

Medical Patient Weight Scale

Model : WM205H

Medical patient weight scale with a maximum capacity of 250 kg. The WM205H comes complete with built in digital height rod and provides automatic calculation of a patient's BMI. Ideal for hospitals, medical clinics and fitness clubs.

FEATURES

WM205H

- professional patient scale with integral digital height rod for determination of weight, height and Body Mass Index (BMI) in a single process
- ideal for hospitals, medical clinics, fitness clubs and retirement homes
- easy-to-read 20 mm LCD with backlight
- auto weight-lock function on stable weight
- weight hold and last height, weight and BMI recall function
- digital height rod with measuring range of 110 cm – 200 cm x 1 mm
- height and BMI automatic calculation and display
- heavy duty load cell platform, with solid robust cast iron construction
- castor wheels at rear for ease of mobility
- AC/DC adaptor and battery operation (6 x AA batteries)
- battery life approximately 75 hours continuous use (alkaline) without backlight
- net weight 15 kg
- Australian Register of Therapeutic Goods (ARTG) # 174843

OPTIONS

- thermal printer providing date, time, weight, height and BMI



Shown with optional printer

SPECIFICATIONS

Model Number	Max Capacity	Graduation	Height Measuring	Platform Size (mm)
WM205H	250 kg	100 g	110 cm - 200 cm x 1 mm	320 (W) x 380 (D) x 90 (H)

Specifications may change without notice